HANDHABUNGSHINWEISE / User Manual

TrayCell

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The fiber-optic ultra-micro cell for UV-Vis analysis

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Contents

1.	Product description	PAGE 5
2.	Product contents	PAGE 6
3.	Product features	PAGE 9
4.	Safety information	PAGE 11
5.	Operation	PAGE 13
6.	Measuring process	PAGE 17
7.	Applications	PAGE 25
8.	Cleaning and maintenance	PAGE 27
9.	Stability of TrayCell materials	PAGE 29
10.	Transportation and storage	PAGE 31
11.	Technical specifications	PAGE 31
12.	Technical specifications: Cap	PAGE 33
13.	FAQs	PAGE 36

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

1. Product description

High tech, tiny footprint, patented operating principle

The TrayCell comprises a fiber-optic measuring cell and a cap with an integrated mirror. A small drop (0.7 μl – 10 μl) of sample is pipetted onto the window and the cap is then placed on top. The precisely defined spacing between the window and the mirror inside the cap ensures that the optical path length is accurate and remains constant. It is

therefore impossible for the path length to change, rendering costly calibrations and readjustments unnecessary. Light is guided through the sample via prisms and fiber-optic waveguides, reflected in the mirror and then guided back out of the TrayCell to the detector via the waveguides.

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- A** 1 × TrayCell with 8.5 mm center height,
type: 105.800-UVS or 105.810-UVS
- B** 1 × cap with 0.2 mm path length [factor 50]*,
Article-No. 665-704-0.2-40
- C** 1 × cap with 1.0 mm path length [factor 10]*,
Article-No. 665-703-1-40
- D** 1 × adapter for 15 mm center height
- E** 1 × adapter for 20 mm center height
- F** 1 × screwdriver for center height adapter
- G** 1 × storage box

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3. Product features

- Strictly defined, precise path lengths ➡ misadjustment impossible
- Path lengths in standard model: 1.0 mm and 0.2 mm
- Additional path lengths: 2.0 mm and 0.1 mm (optional)
- Extremely small sample volumes: 0.7 μl – 10 μl
- Large measuring range from approx. 6 – 8,500 ng/ μl (dsDNA)* and approx. 0.1 – 100 mg/ml (protein)*
- No need to dilute samples
- Cap prevents evaporation of samples
- Easy recovery of samples by pipetting
- Excellent reproducibility
- Quick and easy to clean
- Highly flexible
- Suitable for all current spectrophotometers

Examples for use:

- Determining the purity and concentration of proteins (direct measurement or using colorimetric assay)
- Determining the purity and concentration of DNA/RNA
- Determining labeling efficiency for microarray experiments
- All microvolume, spectrophotometric measurements (0.7 μl – 10 μl) in the UV-Vis range from 190 nm to 1,100 nm

* depending on the spectrophotometer used and the sample type (see example in chapter 6. Measuring range)

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4. Safety information

The TrayCell is solely intended for use in spectrophotometers, e.g. for determining the concentration of analytes in liquids.

The light path of the spectrophotometer is diverted within the TrayCell such that the beam emanating from the light source is able to escape upwards if the mirrored cap is not positioned on the TrayCell.

 **Before each measurement, ensure that the mirrored cap is in place on the TrayCell.**

The TrayCell may not be stored at temperatures below 4°C or above 50°C.

There is a risk of corrosion if aggressive cleaning products and disinfectants are used.

 **Do not use corrosive cleaning agents, aggressive solvents or abrasive polishes.**


The TrayCell consists of quartz glass and metal components, and must therefore be handled with care. After use, return the TrayCell to the storage box and ensure that this is closed. Do not put the TrayCell down in an upright position. If it tips over, the quartz glass optical head may break.

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5. Operation

- The TrayCell is a fiber-optic, microvolume measuring cell and was developed for the UV-Vis analysis of proteins and DNA/RNA. To enable its use in a majority of spectrophotometers, the TrayCell's dimensions are the same as those of a standard cuvette.
- First check the center height required for your spectrophotometer and ensure that the TrayCell is set to the correct center height.
- The TrayCell is supplied with a center height of 8.5 mm. If you need to adjust the center height, please use the relevant adapter. This is supplied with appropriate screws. Screw the adapter into the pre-drilled holes on the TrayCell from below.
- Position the TrayCell in the cuvette holder with the window facing in the direction of the light path. We recommend placing the Hellma logo at the front. Always insert the TrayCell in the same direction.
- Maximum TrayCell performance is achieved when the light beam that passes through the measuring chamber has a diameter below 4 mm.
- The TrayCell must be correctly inserted into the cuvette holder for measurements to be performed. It should then remain in the cuvette holder for the duration of all sample measurements and cleaning processes carried out within a measurement series.
- The TrayCell should not be moved at all (removed and re-inserted) during a series of measurements.

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
Neue Probe pipettieren.


5. Operation

- Please note that the empty TrayCell has an absorbance of around 1. To correct this characteristic absorbance, first carry out a reference measurement without the sample.
- You can then start taking measurements with the sample.
- Please note that the cap must be correctly positioned on the TrayCell before you start measuring!
- **If necessary, the sample can be recovered by pipetting once the measurements are complete.**
- **Finally, clean the window and the mirror inside the cap.**

Insert TrayCell into cuvette holder of the spectrophotometer.

Pipette sample onto window.

Position cap. Start measurement.

Remove cap and recover sample if necessary.

Clean window and cap (TrayCell remains in spectrophotometer).

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6. Measuring range

- By employing a variety of caps, and thus path lengths, the necessity to dilute samples is significantly reduced, which allows the concentration range to be significantly increased as a result.

Example of use: Nucleic acid quantitation

- A solution's ratio of absorbance at 260 nm (A_{260}) is used to determine its nucleic acid content. The following calculation derived from the Beer-Lambert law is used here:

$$\text{Concentration [ng/}\mu\text{l]} = \text{absorbance [260 nm]} \times \text{sample-specific factor} \times \text{virtual dilution factor}$$

- The sample-specific factor indicates the specific absorbance of a material or molecule at a certain wavelength. For example, dsDNA with a concentration of 50 ng/ μ l at 260 nm exhibits a specific absorbance of 1 at a path length of 10 mm. Owing to the varying path lengths possible in the TrayCell, the virtual dilution factor must also be taken into account. This is given as a factor for each cap.

[See listed values for each path length on the next page]

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6. Measuring range

- Absorption measurements at a wavelength of 260 nm on different nucleic acid samples show the mean dynamic range of concentration as a factor of the optical path length and the absorption thresholds of 0.025 – 1.7 as follows:

NUCLEIC ACIDS	SAMPLE-SPECIFIC FACTOR**	2.0 mm cap (factor 5) [ng/μl]*	1.0 mm cap (factor 10) [ng/μl]*	0.2 mm cap (factor 50) [ng/μl]*	0.1 mm cap (factor 100) [ng/μl]*	Total measuring range [ng/μl]*
dsDNA	50	6 – 425	13 – 850	63 – 4250	125 – 8500	6 – 8500
ssDNA	37	5 – 315	9 – 629	46 – 3145	93 – 6290	5 – 6290
ssRNA	40	5 – 340	10 – 680	50 – 3400	100 – 6800	5 – 6800
Oligomers	30	4 – 255	8 – 510	38 – 2550	75 – 5100	4 – 5100

* typical concentration values that can be measured using an average spectrophotometer

** sample-specific factor of nucleic acid samples = concentration [ng/μl] at an absorbance of 1

Sample quantity required	6 – 10 μl	3 – 5 μl	0.7 – 4 μl	0.7 – 3 μl
--------------------------	-----------	----------	------------	------------

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6. Measuring range

Example of use: Protein quantitation

- A solution's absorbance at 280 nm (A_{280}) is often used to determine its protein content. Alternatively, proteins can also be measured through absorbance at 230 nm or in colorimetric assays. The following calculation derived from the Beer-Lambert law is used to perform direct measurements at 280 nm:

$$\text{Concentration [mg/ml]} = \frac{1}{\text{sample specific factor}} \times \text{absorbance [280 nm]} \times \text{dilution factor}$$

- The sample-specific factor indicates the specific absorbance of a protein at a certain wavelength. For example, BSA (bovine serum albumin) with a concentration of 1 mg/ml at 280 nm and a path length of 10 mm exhibits a specific absorbance of 0.64. Owing to the varying path lengths possible in the TrayCell, the virtual dilution factor must also be taken into account. This is given as a factor for each cap.

(See listed values for each path length on the next page)

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6. Measuring range

Absorption measurements at 280 nm on different protein samples show the mean dynamic range of concentration as a factor of the optical path length and the absorption thresholds of 0.025 – 1.0 as follows:

PROTEINS	SAMPLE-SPECIFIC FACTOR**	2.0 mm cap (factor 5) [mg/ml]*	1.0 mm cap (factor 10) [mg/ml]*	0.2 mm cap (factor 50) [mg/ml]*	0.1 mm cap (factor 100) [mg/ml]*	Total measuring range [mg/ml]*
BSA (bovine serum albumin)	0.64	0.2 – 7.8	0.4 – 15.6	2.0 – 78	3.9 – 156	0.2 – 156
IgG (bovine gamma globulin)	1.4	0.1 – 3.6	0.2 – 7.1	0.9 – 36	1.8 – 71	0.1 – 71
Lysozyme	2.64	0.05 – 1.9	0.1 – 3.8	0.5 – 19	0.9 – 38	0.05 – 38

* typical concentration values that can be measured using an average spectrophotometer

** sample-specific factor of protein samples = absorbance (A₂₈₀) at a concentration of 1 mg/ml

Sample quantity required	6 – 10 µl	3 – 5 µl	0,7 – 4 µl	0,7 – 3 µl
--------------------------	-----------	----------	------------	------------

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
7. Applications

TrayCell applications

On our website www.hellma.analytics.com you will find the following link to our application reports, which describe the advantages of the TrayCell for use in spectrophotometers: www.hellma-analytics.com/apps

Through the integrated beam deflection and the use of fibre-optic cables, it is possible to measure the sample directly on the surface of the optical window. The dosing component ("cap") ensures a well-defined optical light path of 1 mm. It comes with openings to prevent pressure build-up when being closed.

Spectral scan

In an example, spectra of DNA in different concentrations were measured. The reference measurement was performed with water. The following parameter settings were used:

Display	Absorbance
Correction	reference
Str [nm]	4
Lamp change at [nm]	320
Meas. mode	step mode
Range [nm]	200 - 350
Delta lambda [nm]	1
Integration time [s]	0.1

Photometric analysis in the ultra-microtitre range
Index: UV_ALL_35_36_v1 (01/2020)_EN

analytikjena
Analytik Jena AG

Photometric analysis in the ultra-microtitre range

General

Photometric measurement of very small sample volumes can be performed with the fibre-optic ultra micro measuring cell and the SPECORD® or SPEKOL® spectrophotometers from Analytik Jena AG. For example, the purity as well as concentration of DNA and RNA can be measured in using just a few microlitres. The analysis of proteins (A280, Bradford, Lowry, Biuret, etc.) can also be performed in the microlitre range. The outstanding photometric properties of the SPEKOL® and SPECORD® series, in combination with the ultra microtitre cell, offer a very broad absorption range.

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The ultra microtitre cell is placed in the SPECOL® / SPECORD® by means of the adjustable holder and is optimally positioned. A centre height of 15 mm or 8.5 mm can be used with the SPECORD® with the relevant holder. The SPEKOL® measures at a centre height of 8.5 mm.

The measuring cell does not need to be removed from the spectrophotometer for filling or cleaning. The guarantee and therefore no variation in positioning of the aperture in the light beam measurement. After the cap is removed, the sample to be measured is directly pipetted onto the surface of the measuring cell. Just 3 to 5 µl of sample is sufficient! When the cap is replaced the sample can be measured. For cleaning the surface solvents as required.

Adjustable holder SPECOL® (left) and SPECORD® (right)

Photometric analysis in the ultra-microtitre range
Index: UV_ALL_35_36_v1 (01/2020)_EN

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8. Cleaning and maintenance

- Always clean the TrayCell with a lint-free cloth or swab.
- Please observe the following when cleaning the TrayCell:

- ➔ Do not immerse the TrayCell in water or cleaning solutions
- ➔ Do not clean the TrayCell in an ultrasonic cleaner
- ➔ Do not use abrasive cleaning agents

- Depending on the sample, use a 60% isopropanol solution, ethanol or ultrapure water to clean the TrayCell. If necessary, the TrayCell can be wiped clean with the solvent used to dissolve the sample.
- Hellmanex III cuvette cleaning concentrate is ideally suited for optimal TrayCell cleaning.
- The TrayCell consists of quartz glass and metal components, and must therefore be handled with care. While the quartz glass components are highly resistant to laboratory cleaning agents, attention must be paid to ensure that no caustic or corrosive cleaning agents are used on metal components.
- After use and cleaning, return the TrayCell to the storage box supplied and ensure that this is closed.

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9. Stability of TrayCell materials

At room temperature, the TrayCell measuring head, including the cap, is resistant to many organic solvents, acids up to $\text{pH} \geq 2$ and alkaline solutions up to $\text{pH} \leq 10$.

The TrayCell measuring head is resistant to the following chemical compounds at room temperature:

- ✓ Acetone (up to 5%)
- ✓ Acetonitrile
- ✓ Benzene
- ✓ Toluene
- ✓ Phenol (up to 1%)
- ✓ Carbon tetrachloride
- ✓ Chloroform
- ✓ Dichloromethane
- ✓ Methanol
- ✓ Ethanol
- ✓ Butanol
- ✓ n-Propanol
- ✓ Isopropyl alcohol
- ✓ Ethers
- ✓ Hexane
- ✓ HEPES
- ✓ MES
- ✓ MOPS
- ✓ Saline buffers, e.g. PBS, Citrate or borate buffer in a pH range from 4 – 10
- ✓ Low (!) concentrations of acids and alkaline solutions

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10. Transport and Storage

- Only transport and store the TrayCell in the storage box supplied.
- Do not store the TrayCell at temperatures below 4°C or above 50°C.

11. Technical specifications

TYPE	105.800-UVS	105.810-UVS
Product No.	105800-A3-V1-46	105810-A3-V1-46
Window material	Quartz SUPRASIL®	Quartz SUPRASIL®
Width x depth	12.5 x 12.5 mm	12.5 x 12.5 mm
Height*	68.5 mm (center height 8.5 mm) 75 mm (center height 15 mm) 80 mm (center height 20 mm)	53 mm (center height 8.5 mm) 59.5 mm (center height 15 mm) 64.5 mm (center height 20 mm)
Volumes	0,7 - 10 µl	0,7 - 10 µl
Path length	0.2 mm or 1 mm (+/- 0.02 mm) 0.1 mm or 2 mm (Optional)	0.2 mm or 1 mm (+/- 0.02 mm) 0.1 mm or 2 mm (Optional)
Max. temperature	50 °C	50 °C
Center height**	8.5 mm, 15 mm or 20 mm* (other center heights available on request)	8.5 mm, 15 mm or 20 mm* (other center heights available on request)
Fiber optics	built in, not exchangeable UV/Vis low solarisation 190 nm – 1,100 nm (52.632 cm ⁻¹ – 9.100 cm ⁻¹)	built in, not exchangeable UV/Vis low solarisation 190 nm – 1,100 nm (52.632 cm ⁻¹ – 9.100 cm ⁻¹)

* The correct height will depend on the cuvette holder design and the spectrophotometer model

The TrayCell should be easily accessible within the cuvette holder of the spectrophotometer and it should be possible to close the cap.

** The center height can be adjusted using the adapter supplied. When ordering, please state the required center height or the brand and model of your spectrophotometer.

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12. Technical specifications: Cap

	STANDARD		OPTIONAL	
Type	665.703	665.704	665.705	665.706
Product No.	665-703-1-40	665-703-0.2-40	665-705-2-40	665-706-0.1-40
Description	Stainless steel cap with integrated mirror for adjusting path lengths			
Path length	1.0 mm	0.2 mm	2.0 mm	0.1 mm
Mirror material	Quartz glass SUPRASIL® with aluminium mirror coating			

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13. FAQs

1. Which TrayCell is right for my device?
2. Where should I pay special attention when inserting the TrayCell?
3. When do I use which cap?
4. How do I clean the TrayCell?
5. In which wavelength range can I take measurements?
6. How do I use the TrayCell in a double beam spectrophotometer?
7. What is the concentration range covered?
8. What is the maximum absorbance I can measure?
9. Can I measure highly-concentrated proteins?
10. Can I measure samples with a low surface tension?
11. My measurement values fluctuate, what can I do?

1. **Which TrayCell is right for my device?** Hellma offers the TrayCell in two different lengths. The 'short' TrayCell 105.810-UVS fits well in devices where the sample compartment is closed with a cap. The 'long' TrayCell 105.800-UVS makes pipetting easier in devices with deep sample compartments. The properties of the two measuring cells are otherwise identical. Spacers can be used in both models to achieve various center heights as necessary. Your TrayCell can therefore be used problem-free in spectrophotometers with different center heights. If you are unsure as to which center height you require, tell us the details of your device and we will recommend the most suitable TrayCell.
2. **Where should I pay special attention when inserting the TrayCell?** As is the case when using cuvettes, you should be careful to ensure that the TrayCell stands straight and stable in the light path. To achieve the best possible reproducibility

of measurement values, we recommend always inserting the TrayCell in the same direction, with the Hellma logo facing the front, and leaving it in the holder between measurements.

3. **When do I use which cap?** The smaller the path length, the higher the concentrations of samples that can be measured within the linearity range of the spectrophotometer or sample testing system. In contrast to the customary cuvettes with their 10 mm path length, the 1 mm path length cap offers a 'dilution factor' of 10, and the cap with a 0.2 mm path length even provides a 'dilution factor' of 50. This factor therefore indicates the extent of dilution that would be necessary in a conventional cuvette [virtual dilution factor].
4. **How do I clean the TrayCell?** The TrayCell consists of quartz glass and metal components, and must therefore be handled with care. While the quartz glass components are highly resistant to laboratory cleaning agents, no caustic or corrosive cleaning agents should be used on metal components. We recommend using lint-free swabs or lint-free lab cloths for cleaning. Samples can be recovered beforehand using a pipette, and the remainder is then wiped away using a swab or lab cloth. If necessary, final cleaning may be performed with a solvent used for the sample. In addition, our Hellmanex III cleaning agent is perfect for thorough TrayCell cleaning.
5. **In which wavelength range can I take measurements?** The TrayCell is characterized by its use of low solarization fiber optics. It permits a measuring range of between 190 and 1,100 nm.
6. **How do I use the TrayCell in a double beam spectrophotometer?** It is not possible to simply use a second TrayCell and follow the usual procedure for taking reference measurements in a double beam spectrophotometer. The differences in the optical characteristics of the fibers are too great to allow any sensible comparisons to be drawn. Baseline correction is sufficient in the ma-

jority of cases, however, making such comparisons unnecessary for standard measurements. If you need to improve the signal-to-noise ratio, we recommend reducing the reference beam to an intensity of around 20%. This can easily be achieved using a simple aperture in the reference beam path.

- 7. What is the concentration range covered?** Depending on the sample being analyzed (double-stranded or single-stranded DNA or RNA, oligomers, etc.) we see nucleic acid concentration ranges from around 6 to 8,500 ng/μl for an absorbance range of 0.025 – 1.7 (see chapter 6. Measuring range). In the case of proteins, measurements should only be taken up to a maximum absorbance of 1, since the scattering above this level would no longer ensure linearity as defined by the Beer-Lambert law. For a protein with an absorbance of 1 at a concentration of 1 mg/ml, the concentration range with a TrayCell would be 0.1 – 100 mg/ml. This concentration range varies depending on the specific absorption coefficients of the individual proteins [see chapter 6. Measuring range].
- 8. What is the maximum absorbance I can measure?** The maximum absorbance that can be measured is limited by the linearity range of the spectrophotometer used. The highest performance spectrophotometers available on the market permit measurements up to a maximum absorbance of 10. Statements in the literature referring to far greater absorbance ratios should be understood as theoretical values. These 'theoretical' absorbance ratios indicate the results that could be expected if we were able to measure undiluted samples in cuvettes with 10 mm path lengths. When using a path length of 0.2 mm (factor 50), the absorbance range of a spectrophotometer with a linear measuring range up to 1.7 would equate to a 'theoretical' absorbance of up to 1.7×50 (factor) = 85 in a 10 mm cuvette.
- 9. Can I measure highly-concentrated proteins?** The ability to measure highly-concentrated samples is determined by the linear absorbance range of the spectrophotometer used. The type of sample being measured also plays a role, as proteins, for example, can only be measured correctly up to an absorbance of 1. By using different path lengths with the TrayCell, even highly-concentrated proteins can be measured.
- 10. Can I measure samples with a low surface tension?** Yes. As the window is slightly indented and measurements are performed in the horizontal plane, the liquid is unable to escape. To achieve the best results, we recommend using the cap with the 0.2 mm path length where possible, or the 0.1 mm path length cap that is available as an optional accessory.
- 11. My measurement values fluctuate, what can I do?** Check that a sufficiently large sample quantity has been pipetted onto the window. The volume should be equivalent to the minimum quantity recommended for a given path length. Some pipettes are not precise enough for transferring small sample volumes. If in doubt, increase the sample volume a little. In the case of very low volumes, we recommend pipetting samples directly onto the mirror inside the cap. Check whether the spectrum is showing a high noise level, that could cause the measurement signal to fluctuate. Provided that the concentration and/or absorbance of the sample lies within the measuring range, an extended integration period will improve the measurements.

NOTIZEN

A series of horizontal dotted lines for taking notes, spanning the width of the page.

Note

A series of 18 horizontal dotted lines for taking notes.

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Hellma GmbH & Co. KG