TRAYCELL 2.0

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We are pleased that you have decided to purchase the TrayCell 2.0 from Hellma. You have acquired a high-precision product with simple handling and highest Hellma quality. This quick start guide provides you with the most important information on the correct handling of the product. For more detailed information please download the complete manual on our website. ► www.hellma.com/traycell-en
If you have any further questions, you can also contact us personally. ► sales.analytics@hellma.com

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Please find the picture with delivery content on the other page.

Optional Accessories

- 1 x cap with 0.1 mm path length (factor 100)*,
Article-No. 665-1016-0.1-40
- 1 x cap with 2.0 mm path length (factor 5)*,
Article-No. 665-1016-2-40

* based on 10 mm path length

SAFETY INFORMATION

The TrayCell is solely intended for use in spectrophotometers, e.g. for determining the concentration of analytes in liquids.

The light path of the spectrophotometer is diverted within the TrayCell such that the beam emanating from the light source is able to escape upwards if the mirrored cap is not positioned on the TrayCell.

 **Before each measurement, ensure that the mirrored cap is in place on the TrayCell.**

The TrayCell must not be stored at temperatures below 4°C or above 50°C.

There is a risk of corrosion if aggressive cleaning products and disinfectants are used.

 **Do not use corrosive cleaning agents, aggressive solvents or abrasive polishes.**

The TrayCell consists of quartz glass and metal components, and must therefore be handled with care. After use, return the TrayCell to the storage box and ensure that this is closed. The TrayCell can easily tip over in an upright position, which can result in damage. You should therefore always place the TrayCell on a lint-free surface or in the storage box.


OPERATION

- The TrayCell is a fiber-optic, microvolume measuring cell and was developed for the UV-Vis analysis of proteins and DNA/RNA. To enable its use in a majority of spectrophotometers, the TrayCell's dimensions are the same as those of a standard cuvette.
- First check the center height required for your spectrophotometer and ensure that the TrayCell is set to the correct center height.
- Position the TrayCell in the cuvette holder with the window facing in the direction of the light path. We recommend placing the Hellma logo at the front. Always insert the TrayCell in the same direction.
- Maximum TrayCell performance is achieved when the light beam that passes through the measuring chamber has a diameter below 4 mm.
- The TrayCell must be correctly inserted into the cuvette holder for measurements to be performed. Please ensure that the TrayCell is firmly seated in the cuvette holder and does not wobble! It should then remain in the cuvette holder for the duration of all sample measurements and cleaning processes carried out within a measurement series.
- The TrayCell should not be moved at all (removed and re-inserted) during a series of measurements.


The TrayCell is supplied with a center height of 8.5 mm. If you need to adjust the center height, please use the relevant adapter. This is supplied with appropriate screws. Screw the adapter into the pre-drilled holes on the TrayCell from below.

- Please note that the empty TrayCell has an absorbance of around 1. To correct this characteristic absorbance, first carry out a reference measurement with the solvent in which the sample was dissolved.
- You can then start taking measurements with the sample.
- Please note that the cap must be correctly positioned on the TrayCell before you start measuring!
- If necessary, the sample can be recovered by pipetting once the measurements are complete.
- Finally, clean the window and the mirror inside the cap.


A Pipette sample onto window.

B Position cap.
Start measurement.

Attention: Place the cap in such a way that the pins are in the notches of the cap.

C Remove cap and recover sample if necessary.

D Clean window and cap (TrayCell remains in spectrophotometer).

E Pipette new sample onto window.

MEASURING RANGE

- By employing a variety of caps (path lengths), and thus virtual dilution factors, the necessity to dilute samples is significantly reduced, which allows the concentration range to be significantly increased as a result.

Example of use: Nucleic acid quantitation

- A solution's ratio of absorbance at 260 nm (A260) is used to determine its nucleic acid content. The following calculation derived from the Beer-Lambert law is used here:

$$\text{Concentration [ng/}\mu\text{l]} = \text{absorbance (260 nm)} \times \text{sample-specific factor} \times \text{virtual dilution factor}$$

Example of use: Protein quantitation

- A solution's absorbance at 280 nm (A280) is often used to determine its protein content. Alternatively, proteins can also be measured through absorbance at 230 nm or in colorimetric assays. The following calculation derived from the Beer-Lambert law is used to perform direct measurements at 280 nm:

$$\text{Concentration [mg/ml]} = \frac{1}{\text{sample specific factor}} \times \text{absorbance (280 nm)} \times \text{virtual dilution factor}$$

On our website you can download a helpful Excel tool, with which you can calculate the concentrations of nucleic acids and proteins:

www.hellma.com/traycell-en

STABILITY OF TRAYCELL MATERIALS

At room temperature, the TrayCell measuring head, including the cap, is resistant to many organic solvents, acids up to $\text{pH} \geq 2$ and alkaline solutions up to $\text{pH} \leq 10$. The TrayCell measuring head is resistant to the following chemical compounds at room temperature:

- ✓ Acetone (up to 5%)
- ✓ Acetonitrile
- ✓ Benzene
- ✓ Toluene
- ✓ Phenol (up to 1%)
- ✓ Carbon tetrachloride
- ✓ Chloroform
- ✓ Dichloromethane
- ✓ Methanol
- ✓ Ethanol
- ✓ Butanol
- ✓ n-Propanol
- ✓ Isopropyl alcohol
- ✓ Ethers
- ✓ Hexane
- ✓ HEPES
- ✓ MES
- ✓ MOPS
- ✓ Saline buffers, e.g. PBS,
Citrate or borate buffer in a pH range from 4 – 10
- ✓ Low (!) concentrations of acids and alkaline solutions

WE ARE PLEASED
TO OFFER HELP AND ADVICE.

Hellma GmbH & Co. KG

Klosterrunsstraße 5

79379 Müllheim

Germany

phone + 49 7631 182 0

info.de@hellma.com

www.hellma.com

Find all international contacts on:
www.hellma.com/contacts